

# COVID-19 and Implications in Haemostasis and Viscoelastic Testing

## What is COVID-19?

Coronavirus Disease 2019 (COVID-19) is an infection caused by the novel coronavirus, SARS-CoV-2, and is characterised by flulike symptoms, such as fever and cough. Some patients may develop dyspnea, tachypnea and/or pneumonia with disturbed gas exchange.<sup>1</sup> In severe cases, the disease may lead to acute respiratory distress syndrome (ARDS), sepsis and multiorgan failure.<sup>2</sup>

**COVID-19 severity can be classified as mild, moderate, severe or critical.<sup>3-5</sup>**

Classification	Symptoms
Asymptomatic	None
Mild	Mild clinical symptoms and no signs of pneumonia
Moderate	Fever and respiratory symptoms
Severe	Patients with any of the following: respiratory distress with respiratory rate $\geq 30$ breaths/min, $SpO_2 \leq 93\%$ at rest, $PaO_2/FiO_2 \leq 300$ mm Hg (mmHg = 0.133 kPa)
Critical	Patients with any of the following: respiratory failure requiring mechanical ventilation, shock or other organ failure requiring admission to Intensive Care Unit (ICU)

## Comorbidities and COVID-19

Patients with certain comorbidities may be at higher risk for severe SARS-CoV-2 infection. The following disorders have been shown to be associated with increased COVID-19 severity.<sup>6</sup>

- Hypertension
- Cerebrovascular disease
- Diabetes
- Chronic obstructive pulmonary disease
- Cardiovascular disease


The information in this document is accurate as of posting time. However, as the situation surrounding COVID-19 continues to evolve, it's possible that some data will have changed since publication. While IL is trying to keep COVID-19 content as up to date as possible, we also encourage readers to stay informed on news and recommendations by referencing guidelines published by the ISTH, ASH, CDC, and WHO and consulting with local healthcare authorities.

# COVID-19 and Haemostasis

Early reports indicate that patients with COVID-19 may develop abnormalities in haemostasis laboratory testing.<sup>5,7,8</sup> In a recent meta-analysis, elevations of D-dimer were reported in 20.4% of COVID-19 patients.<sup>9</sup> These abnormalities have subsequently been correlated with disease severity and outcomes including mortality.<sup>10–23</sup> Moreover, thrombotic complications including venous thromboembolism (VTE) have been reported in 16–42.7% of patients with COVID-19.<sup>24–26</sup> The pathogenesis of the coagulopathy associated with COVID-19 is an area of active investigation. Preliminary hypotheses suggest hyperinflammation and associated immunothrombosis may be responsible.<sup>27–29</sup>

## Role of Haemostasis Testing in Coagulopathy and COVID-19

Haemostasis testing provides important information for the management of COVID-19, especially in severe cases. The International Society on Thrombosis and Haemostasis (ISTH) “interim guidance on recognition and management of coagulopathy in COVID-19” recommends the use of haemostasis assays, D-dimer, prothrombin time (PT), platelet count and fibrinogen to assess prognosis in COVID-19 patients who require admission.<sup>30</sup>



Algorithm for the management of coagulopathy in COVID-19 based on simple laboratory markers.

\*The list of markers is given in decreasing order of importance.

\*\*Performing fibrinogen assays may not be feasible in many laboratories but monitoring the levels can be helpful after patient admission.

\*\*\*Although a specific cut-off cannot be defined, a three- to four-fold increase in D-dimer values may be considered significant. Any one of the values in this table may be considered significant.

**Figure 1** Adapted from: Thachil J, Tang N, Gando S, *et al.* ISTH interim guidance on recognition and management of coagulopathy in COVID-19. *J Thromb Haemost.* 2020;18(5):1023–6. doi:10.1111/jth.14810.

# Impact of COVID-19 on Haemostasis Markers

Publications have reported the following changes in haemostasis assays in patients with COVID-19:

Assay	COVID-19 impact
APTT	Normal <sup>7, 9, 11, 16-20, 26, 31-36</sup>
PT	Normal to ↑ <sup>7, 9, 11, 16-20, 26, 31-36</sup>
Thrombin Time	Normal to ↑ <sup>11, 33</sup>
Fibrinogen	↑↑ <sup>11, 17, 18, 20, 26, 32, 34-36</sup>
Platelet Count	↓ <sup>37, 38</sup>
D-Dimer	↑↑↑ <sup>7, 9, 11, 13-20, 26, 31-36, 39-41</sup>
FDP	↑↑ <sup>11, 18, 32</sup>
Antithrombin	Normal to ↓ <sup>11, 26, 32, 34-36</sup>
Protein C	Normal to ↑ <sup>36</sup>
Free Protein S Antigen	Normal to ↓ <sup>11, 32, 34-36</sup>
Factor V Activity	Normal <sup>26</sup>
Factor VIII Activity	↑↑ <sup>26, 36, 42</sup>
VWF:Act, VWF:Ag or VWF:RCo <sup>15</sup>	Normal to ↑↑ <sup>26, 36, 42</sup>
Lupus Anticoagulant	Normal to Present <sup>26, 43-46</sup>
Anti-Cardiolipin IgA/IgM/IgG	Normal to ↑ <sup>43, 47</sup>
Anti-β <sub>2</sub> Glycoprotein I IgA/IgM/IgG	Normal to ↑ <sup>43, 47</sup>
ADAMTS13 Activity	Normal to ↓ <sup>11, 32, 34-36, 42</sup>

APTT = activated partial thromboplastin time

FDP = fibrin degradation products

VWF:Act: von Willebrand Factor Activity

VWF:Ag: von Willebrand Factor Antigen

VWF:RCo: von Willebrand Factor Ristocetin Cofactor Activity

## Viscoelastic Testing

Early reports also suggest hypercoagulability in patients with COVID-19, defined by short clot formation times and increased clot firmness from viscoelastic testing.<sup>34–36, 48, 49</sup> In addition, there has not been evidence of secondary hyperfibrinolysis, supporting severe hypercoagulability, rather than a consumptive coagulopathy as part of the pathogenesis of SARS-CoV-2 infection.<sup>48</sup> These preliminary findings indicate that viscoelastic testing may have a role in the rapid identification of patients with severe COVID-19. Moreover, viscoelastic testing may have clinical utility in assessing a patient's response to anticoagulant therapy.

## Thrombotic Complications of COVID-19 and Anticoagulation Therapy


Several studies have demonstrated a high incidence (16–43%) of thrombotic complications in COVID-19 patients<sup>24</sup>, including:

- VTE
- Thrombosis of extracorporeal circuits (e.g., extracorporeal membrane oxygenation, continuous renal replacement therapy)
- Central venous catheter-associated thrombosis

These thrombotic complications occurred, despite prophylactic or therapeutic anticoagulation, and were more likely to occur in critically ill COVID-19 patients. Despite this apparent failure of anticoagulation therapy in some patients, early studies indicate that the use of low molecular weight heparin (LMWH) is associated with reduced mortality in COVID-19 patients with coagulopathy.<sup>50</sup>

Current recommendations advise treating all hospitalised COVID-19 patients with standard-dose LMWH, unfractionated heparin (UFH) or fondaparinux for VTE prophylaxis, unless contraindicated (e.g., bleeding, thrombocytopenia).<sup>1, 3, 51–53</sup> The use of intermediate dose LMWH should be considered for patients with multiple risk factors for VTE (e.g., immobilisation, obesity, history of VTE).<sup>1</sup>

Elevated FVIII activity and the presence of lupus anticoagulants make using APTT for heparin monitoring challenging. Thus, use of an anti-Xa assay to monitor UFH therapy in these patients is recommended.<sup>3</sup> Heparin-induced thrombocytopenia (HIT) is a risk in all patients on heparin therapy. Regular monitoring and 4T scores should be assessed and if HIT is present, patients should be converted to a non-heparin anticoagulant therapy (e.g., danaparoid, argatroban, bivalirudin).<sup>54, 55</sup>



**Figure 2** Adapted from: Bikdeli B, Madhavan MV, Jimenez D, *et al.* COVID-19 and thrombotic or thromboembolic disease: Implications for prevention, antithrombotic therapy, and follow-up [published online ahead of print, Apr 15, 2020]. *J Am Coll Cardiol.* 2020. doi:10.1016/j.jacc.2020.04.031.

## References

- Marietta M, Ageno W, Artoni A, *et al.* COVID-19 and haemostasis: a position paper from Italian Society on Thrombosis and Haemostasis (SISST) [published online ahead of print, April 8, 2020]. *Blood Transfus.* 2020. doi:10.2450/2020.0083-20.
- World Health Organization. Clinical management of severe acute respiratory infection (SARI) when COVID-19 disease is suspected: interim guidance. March 2020. WHO/2019-nCoV/clinical/2020.5.
- Zhai Z, Li C, Chen Y, *et al.* Prevention and treatment of venous thromboembolism associated with coronavirus disease 2019 infection: A consensus statement before guidelines [published online ahead of print, Apr 21, 2020]. *Thromb Haemost.* 2020. doi:10.1055/s-0040-1710019.
- Yang X, Yu Y, Xu J, *et al.* Clinical course and outcomes of critically ill patients with SARS-CoV-2 pneumonia in Wuhan, China: a single-centered, retrospective, observational study [published online ahead of print, 2020 Feb 24] [published correction appears in *Lancet Respir Med.* Apr 8, 2020 (4):e26]. *Lancet Respir Med.* 2020. doi:10.1016/S2213-2600(20)30079-5.
- Guan W, *et al.* Clinical characteristics of coronavirus disease 2019 in China. *N Engl J Med.* 2020;382:1708-20.
- Wang B, Li R, Lu Z, Huang Y. Does comorbidity increase the risk of patients with COVID-19: evidence from meta-analysis. *Aging (Albany, NY).* 2020;12(7):6049-57. doi:10.18632/aging.103000.
- Wang D, Hu B, Hu C. Clinical characteristics of 138 hospitalised patients with 2019 novel coronavirus-infected pneumonia in Wuhan, China. *JAMA.* 2020;323(11):1061-9.
- Huang C, Wang Y, Li X, *et al.* Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China. *Lancet* 2020;395(10223):497-506.
- Henry BM, de Oliveira MHS, Benoit S, *et al.* Hematologic, biochemical and immune biomarker abnormalities associated with severe illness and mortality in coronavirus disease 2019 (COVID-19): a meta-analysis. [published online ahead of print, April 10, 2020]. *Clin Chem Lab Med.* 2020. doi:10.1515/ccm-2020-0369.
- Zhou F, Yu T, Du R, *et al.* Clinical course and risk factors for mortality of adult inpatients with COVID-19 in Wuhan, China: a retrospective cohort study. *Lancet.* 2020;395(10229):1054-62.
- Han H, Yang L, Liu R, *et al.* Prominent changes in blood coagulation of patients with SARS-CoV-2 infection. *Clin Chem Lab Med.* 2020;58(7):1116-20.
- Wang D, Yin Y, Hu C, *et al.* Clinical course and outcome of 107 patients infected with the novel coronavirus, SARS-CoV-2, discharged from two hospitals in Wuhan, China. *Crit Care.* 2020;24(1):188.
- Yao Q, Wang P, Wang X, *et al.* Retrospective study of risk factors for severe SARS-CoV-2 infections in hospitalized adult patients [published online ahead of print, May 29, 2020]. *Pol Arch Intern Med.* 2020. doi:10.20452/pamw.15312.
- Zhang J, Liu P, Wang M, *et al.* The clinical data from 19 critically ill patients with coronavirus disease 2019: a single-centered, retrospective, observational study [published online ahead of print, April 21, 2020]. *Z Gesundh Wiss.* 2020. doi:10.1007/s10389-020-01291-2.
- Zhang JJ, Dong X, Cao YY, *et al.* Clinical characteristics of 140 patients infected with SARS-CoV-2 in Wuhan, China [published online ahead of print, Feb 19, 2020]. *Allergy.* 2020. doi:10.1111/all.14238.
- Li Y, Hu Y, Yu J, Ma T. Retrospective analysis of laboratory testing in 54 patients with severe- or critical-type 2019 novel coronavirus pneumonia. [published online ahead of print, April 27, 2020]. *Lab Invest.* 2020. doi:10.1038/s41374-020-0431-6.
- Zhang L, Yan X, Fan Q, *et al.* D-dimer levels on admission to predict in-hospital mortality in patients with COVID-19. *J Thromb Haemost.* 2020;18:1324-9.
- Zou Y, Guo H, Zhang Y, *et al.* Analysis of coagulation parameters in patients with COVID-19 in Shanghai, China. [published online ahead of print, April 30, 2020]. *Biosci Trends.* 2020. doi:10.5582/bst.2020.03086.
- Wu C, Chen X, Cai Y, *et al.* Risk factors associated with acute respiratory distress syndrome and death in patients with coronavirus disease 2019 pneumonia in Wuhan, China. [published online ahead of print, March 13, 2020]. *JAMA Intern Med.* 2020. doi:10.1001/jamainternmed.2020.0994.
- Fogarty H, Townsend L, Ni Cheallaigh C, *et al.* COVID-19 coagulopathy in caucasian patients. [published online ahead of print, 2020 Apr 24]. *Br J Haematol.* 2020. doi:10.1111/bjh.16749.
- Cummings MJ, Baldwin MR, Abrams D, *et al.* Epidemiology, clinical course, and outcomes of critically ill adults with COVID-19 in New York City: a prospective cohort study. [published online ahead of print, May 19, 2020]. *Lancet.* 2020. doi:10.1101/2020.04.15.20067157.
- Petrilli CM, Jones SA, Yang J, *et al.* Factors associated with hospital admission and critical illness among 5279 people with Coronavirus disease 2019 in New York City: Prospective Cohort Study. *BMJ.* 2020. doi:10.1136/bmj.m1966.
- Cecconi M, Piovani D, Brunetta E. Early Predictors of Clinical Deterioration in a Cohort of 239 Patients Hospitalised for COVID-19 Infection in Lombardy, Italy. *J Clin Med.* 2020. doi:10.3390/jcm9051548.
- Middelorp S, Coppens M, van Haaps TF, *et al.* Incidence of venous thromboembolism in hospitalised patients with COVID-19 [published online ahead of print, May 5, 2020]. *J Thromb Haemost.* 2020. doi:10.1111/jth.14888.
- Klok FA, Kruip MJHA, van der Meer NJM, *et al.* Incidence of thrombotic complications in critically ill ICU patients with COVID-19. *Thromb Res.* 2020;191:145-7. doi:10.1016/j.thromres.2020.04.013.
- Helms J, Tacquard C, Severac F, *et al.* High risk of thrombosis in patients with severe SARS-CoV-2 infection: a multicenter prospective cohort study [published online ahead of print, May 4, 2020]. *Intensive Care Med.* 2020;1-10. doi:10.1007/s00134-020-06062-x.
- Henry BM, Vikse J, Benoit S, *et al.* Hyperinflammation and derangement of renin-angiotensin-aldosterone system in COVID-19: A novel hypothesis for clinically suspected hypercoagulopathy and microvascular immunothrombosis. *Clin Chim Acta.* 2020;507:167-73. doi:10.1016/j.cca.2020.04.027.
- Thachil J, Agarwal S. Understanding the COVID-19 coagulopathy spectrum. *Anaesthesia.* 2020. doi:10.1111/anae.15141.
- Iba T, Levy JH, Levi M, *et al.* Coagulopathy of coronavirus disease 2019 [published online ahead of print, May 27, 2020]. *Crit Care Med.* 2020. doi:10.1097/CCM.0000000000004458.
- Thachil J, Tang N, Gando S, *et al.* ISTH interim guidance on recognition and management of coagulopathy in COVID-19. *J Thromb Haemost.* 2020;18(5):1023-6. doi:10.1111/jth.14810.
- Xiong M, Liang X, Wei YD. Changes in blood coagulation in patients with severe coronavirus disease 2019 (COVID-19): a meta-analysis [published online ahead of print, Apr 18, 2020]. *Br J Haematol.* 2020. doi:10.1111/bjh.16725.
- Tang N, Li D, Wang X, Sun Z. Abnormal coagulation parameters are associated with poor prognosis in patients with novel coronavirus pneumonia. *J Thromb Haemost.* 2020;18(4):844-7. doi:10.1111/jth.14768.
- Wang K, Zuo P, Liu Y, *et al.* Clinical and laboratory predictors of in-hospital mortality in patients with COVID-19: a cohort study in Wuhan, China [published online ahead of print, May 3, 2020]. *Clin Infect Dis.* 2020. doi:10.1093/cid/ciaa538.
- Spiezia L, Boscolo A, Poletto F, *et al.* COVID-19-related severe hypercoagulability in patients admitted to intensive care unit for acute respiratory failure. *Thromb Haemost.* 2020;120(06):998-1000.
- Ranucci M, Ballotta A, Di Dedda U, *et al.* The procoagulant pattern of patients with COVID-19 acute respiratory distress syndrome [published online ahead of print, Apr 17, 2020]. *J Thromb Haemost.* 2020. doi:10.1111/jth.14854.
- Panigada M, Bottino N, Tagliabue P, *et al.* Hypercoagulability of COVID-19 patients in intensive care unit. A Report of Thromboelastography Findings and other Parameters of Haemostasis [published online ahead of print, Apr 17, 2020]. *J Thromb Haemost.* 2020. doi:10.1111/jth.14850.
- Yang X, Yang Q, Wang Y, *et al.* Thrombocytopenia and its association with mortality in patients with COVID-19. *J Thromb Haemost.* 2020;18: 1469-1472. doi.org/10.1111/jth.14848.
- Lippi G, Plebani M, Henry BM. Thrombocytopenia is associated with severe coronavirus disease 2019 (COVID-19) infections: A meta-analysis. *Clin Chim Acta.* 2020;506:145-8. doi:10.1016/j.cca.2020.03.022.
- Terpos E, Ntanasis-Stathopoulos I, Elalamy I, *et al.* Haematological findings and complications of COVID-19 [published online ahead of print, Apr 13, 2020]. *Am J Hematol.* 2020. doi:10.1002/ajh.25829.
- Gao Y, Li T, Han M, *et al.* Diagnostic utility of clinical laboratory data determinations for patients with the severe COVID-19 [published online ahead of print, Mar 17, 2020]. *J Med Virol.* 2020. doi:10.1002/jmv.25770.
- Zhu J, Zhong Z, Ji P, *et al.* Clinicopathological characteristics of 8697 patients with COVID-19 in China: a meta-analysis [published correction appears in *Fam Med Community Health.* May 8, 2020]. *Fam Med Community Health.* 2020. doi:10.1136/fmch-2020-000406.
- Zheng XL. ADAMTS13 and von Willebrand factor in thrombotic thrombocytopenic purpura. *Annu Rev Med.* 2015;66:211-25. doi:10.1146/annurev-med-061813-013241.
- Harzallah I, Debliquis A, Drénou B. Lupus anticoagulant is frequent in patients with COVID-19 [published online ahead of print, Apr 23, 2020]. *J Thromb Haemost.* 2020. doi:10.1111/jth.14867.
- Tang N. Response to "Lupus anticoagulant is frequent in patients with COVID-19" [published online ahead of print, May 7, 2020]. *J Thromb Haemost.* 2020. doi:10.1111/jth.14890.
- Escher R, Breakey N, Lämmle B. Severe COVID-19 infection associated with endothelial activation. *Thromb Res.* 2020;190:62. doi:10.1016/j.thromres.2020.04.014.
- Bowles L, Platon S, Yartey N, *et al.* Lupus Anticoagulant and Abnormal Coagulation Tests in Patients with COVID-19 [published online ahead of print, May 5, 2020]. *N Engl J Med.* 2020. doi:10.1056/NEJMc2013656.
- Zhang Y, Xiao M, Zhang S, *et al.* Coagulopathy and Antiphospholipid Antibodies in Patients with COVID-19. *N Engl J Med.* 2020;382(17):e38. doi:10.1056/NEJMc2007575.
- Pavoni V, Giansello L, Pazzi M, *et al.* Evaluation of coagulation function by rotation thromboelastometry in critically ill patients with severe COVID-19 pneumonia [published online ahead of print, May 11, 2020]. *J Thromb Thrombolysis.* 2020. doi:10.1007/s11239-020-02130-7.
- Raval JS, Burnett AE, Rollins-Raval MA, *et al.* Viscoelastic testing in COVID-19: a possible screening tool for severe disease? [published online ahead of print, May 6, 2020]. *Transfusion.* 2020. doi:10.1111/trf.15847.
- Tang N, Bai H, Chen X, *et al.* Anticoagulant treatment is associated with decreased mortality in severe coronavirus disease 2019 patients with coagulopathy. *J Thromb Haemost.* 2020;18(5):1094-9.
- Barnes GD, Burnett A, Allen A, *et al.* Thromboembolism and anticoagulant therapy during the COVID-19 pandemic: interim clinical guidance from the anticoagulation forum [published online ahead of print, May 21, 2020]. *J Thromb Thrombolysis.* 2020. doi:10.1007/s11239-020-02138-z.
- Song JC, Wang G, Zhang W, *et al.* Chinese expert consensus on diagnosis and treatment of coagulation dysfunction in COVID-19. *Mil Med Res.* 2020;7(1):19.
- Spyropoulos AC, Levy JH, Ageno W, *et al.* Scientific and Standardisation Committee communication: clinical guidance on the diagnosis, prevention and treatment of venous thromboembolism in hospitalized patients with COVID-19. [published online ahead of print, May 27, 2020]. *J Thromb Haemost.* 2020. doi:10.1111/jth.14929.
- Liu X, Zhang X, Xiao Y, *et al.* Heparin-induced thrombocytopenia is associated with a high risk of mortality in critical COVID-19 patients receiving heparin-involved treatment. *Lancet.* 2020. doi:10.1101/2020.04.23.20076851.
- Cuker A, Arepally GM, Chong BH, *et al.* American Society of Haematology 2018 guidelines for management of venous thromboembolism: heparin-induced thrombocytopenia. *Blood Adv.* 2018;2(22):3360-92. doi:10.1182/bloodadvances.2018024489.