

## Technical Specifications

**Test menu**..... Clotting, chromogenic and immunological assays

**Continuous sample and reagent loading**..... YES

**Continuous operation**..... YES

**Sample and reagent transport system**..... Racks

**Samples onboard**..... 80 (10/sample rack)

### Assay-specific

**HIL sample check**..... YES

**Tube-fill height detection**... YES

**Sample aspiration clog check**..... YES

**Closed-Tube Sampling**..... YES (cap-piercing)

**Reagents onboard**..... 40 (refrigerated)

**Reagents barcode reader**... YES (integrated)

**Barcoded reagents**..... YES (barcode includes lot, exp. date and vial size)

**Cuvettes onboard**..... 800

**Applications onboard**..... 500 (250 user-definable)

**Tests/sample**..... 30

### Throughput

PT..... up to 240 tests/hr (240 samples/hr)

APTT..... up to 180 tests/hr (180 samples/hr)

PT/APTT..... up to 180 tests/hr (90 samples/hr)

**Sample predilution**..... YES

**Calibration curve predilution**..... YES

**Factor parallelism**..... YES

**STAT capability**..... YES (at any time, in any position)

**Rerun and reflex testing**... YES (configurable)

**Reaction-curves display**... YES

**Quality Control program**... YES (with configurable multi-rules)

**Results autovalidation**..... YES

**Patient sample results database**..... 20,000 samples (configurable)

**Security system**..... YES (configurable)

**Auto log-off**..... YES

**Events log system**..... YES

**Audit-trail report**..... YES

**Electronic signature**..... YES

**Bidirectional interface**..... YES (with host query)

**PC**..... External (dedicated)

**Operating system**..... Windows® 7-based

**Monitor**..... 17-inch color LCD touchscreen (external)

**Keyboard**..... YES (external)

**Mouse**..... YES (external)

**Printer**..... YES (external)

**User interface**..... Windows 7-based

### Unit dimensions

**(w x d x h)**..... 110 x 82 x 73 cm (43 x 32 x 29 in)

**Weight (analyzer)**..... 147 kg (324 lbs)